

Center for Skeletal Research
MGH Endocrine Unit

Bone Cells Core

Director: Paola Divieti Pajevic, M.D., Ph.D.

Elizabeth Petit, Technician

Email: Bonecellcore@partners.org

Website: <https://csr.mgh.harvard.edu/skeletal-cells>

Bone Marrow Cultures

Flush bone marrow from femur/ tibia (cut both end and flush the BM with 1 ml of α -MEM ice cold)

- Collect cells and count cells (take 100-200 μ l of cell suspension, add equal volume of 2% acetic acid to lyse red blood cells. Let sit for 10 min @ room temp, count-remember that cells were diluted 1:2 !) and resuspend the cells in α MEM 10%FBS 1%PS

FOR Osteoclasts formation

- Resuspend the cells to a final density of 1.5×10^6 cells/ml Plate cells at 1.5×10^6 cells/ml in 24-well plates (0.5 ml/well). Use only the 8 central wells and add PBS or H₂O to the outside wells. Cells are plated in α -MEM+10% FBS+1% PS and Dexamethasone 100nM (protocol from Takahasi et al. Endocrinology 1988)
- Add treatment : vit D (10nM) or PTH (10nM) to induce OCs differentiation
- Every other day change $\frac{1}{2}$ volume (250 μ l) and replace with fresh medium containing 2X Dexamethasone and treatment
- Ocs will form within 2 weeks

FOR CFU-s experiment

- Resuspend the cells to a final density of 1×10^6 cells/ml and plate into a 6-well plate (2ml/well), Add mineralization medium (50 μ g/ml of L-Ascorbic Acid and 10mM β -glycerophosphate) for CFU-osteoblasts
- Change medium 3 times/week
- Stain for ALK Phosph and count colonies bigger than 2 mm (Cfu-Ob)
- Stain with crystal violet for total cells count (CFU-f)