

Protocol for IHC of phosphoERK in growth plate
(with tyramide amplification, Perkin Elmer, Cat. No. NEL700)

- Controls:
1. no 1st Ab, or 1:100 nonspecific IgG
 2. no amplification

Preparation of reagents:

1. 90% methanol, 10% DMSO (for 200ml need 180ml methanol, 20ml DMSO)
2. 1l 1x PBS
3. 0.3% H₂O₂ (for 200ml need 2ml 30% H₂O₂ in 200ml methanol)
4. Biotinyl Tyramide: add DMSO (0.3 ml [NEL700A] or 1.2 ml [NEL700]), store at 4C (stable 6 months)
5. Biotinyl Tyramide Working Solution: dilute stock 1:50 with 1x Amplification Diluent, 100ul/slide
6. Wash Buffer (TNT; 0.1M Tris, 0.15M NaCl, 0.05% Tween, pH 7.5): for 3l need 300ml 1M Tris, pH 7.5, 90ml 5M NaCl, 1.5ml Tween, 2610ml water
7. Blocking Buffer (TNB; 0.1M Tris, 0.15M NaCl, 0.5% Blocking Reagent, pH 7.5): for 10ml need 1ml 1M Tris, 0.3ml 5M NaCl, 0.05g Blocking Reagent, heat to 60C to completely dissolve Blocking Reagent, can be stored at -20C in aliquots, at RT no longer than 24h
8. SA-HRP (1:100 in TNB buffer)
9. DAB kit for visualization

Procedural notes:

1. need humidified chamber (slide box with damp paper towels) for ALL incubation steps
2. staining tubs for incubations need volumes of 200ml

Protocol:

1. heat paraffin slides for 30min at 55C to melt wax
2. deparaffinize slides in 2 changes of xylene x 5min
3. rehydrate slides in series of graded alcohol 100%, 80%, 70% x 5 min
4. incubate slides in 90% methanol, 10% DMSO for 20min at RT
5. wash with PBS
6. incubate with 0.3% H₂O₂ for 20min *in methanol*
7. wash with PBS
8. block slides with 100ul TNB for 30min in humidified chamber at RT, drain off
9. apply 100ul of 1st Ab (1:500 anti-phosphoERK in TNB buffer, Cell Signaling, Cat.No. 4377S) for *overnight at 4°C* *at 4°C*
10. wash 3x5min with TNT buffer at RT with agitation
11. apply 100ul 2nd biotinylated Ab (1:300 anti-rabbit biotinylated IgG in TNB, Vector, Cat.No. BA1000) for 1h at RT
12. wash 3x5min with TNT buffer
13. apply 100ul SA-HRP (diluted 1:100 in TNB) for 30min at RT
14. wash 3x5min with TNT buffer
15. apply 100ul Biotinyl Tyramide Working Solution for *5* min at RT

16. wash 3x5min with TNT buffer
17. add 100ul of SA-HRP and incubate 30min at RT
18. wash 3x5min with TNT buffer
19. visualize with chromogenic substrate 5min in the DARK
 - 5ml H₂O
 - 4 drops DAB
 - 2 drops buffer
 - 2 drops H₂O₂
 - 2 drops NiCl (optional: will give grey/black stain as compared to brown with DAB only) ^{YES}*look for signal while developing*
20. wash with distilled water until desired signal intensity achieved (1-3min)
21. counterstain with hematoxylin (5sec; Histo-core sequence) or toluidine blue (see separate protocol)